



**Web  
Scientific**

# Benchtop Cryo-Station

## Keep Liquid N<sub>2</sub> On Your Bench For Hours!

Introducing the latest in Cryogenic sample handling - the **Benchtop Cryo-Station** allows the rapid, efficient freezing of samples in a robust breakage-free container, specifically designed to retain your liquid N<sub>2</sub> for as long as possible. Combine with the Cryo Dipper tube holder or Cryo Net for added convenience and versatility.



### Benchtop Cryo-Station

- Super insulating triple density expanded polystyrene
- Low thermal mass of EPS minimises liquid N<sub>2</sub> boil-off
- Lightweight construction
- No glass liner to break
- 2 inch thick lid eliminates condensation



### Cryo Dipper

- Quickly load tubes for freezing
- Easily retrieve tubes
- Holds cryogenic vials, 0.5 - 2.0ml
- Movable weights prevent floating
- Durable PVC construction



### Cryo Net

- For quick sample freezing
- Remove all tubes at once
- Durable nylon mesh
- Open design for unimpeded liquid N<sub>2</sub> access

### Ordering Information

Cat. No.	Description	Pack Size	List Price
CRST-1000	Benchtop Cryo-Station	1	£98
CDRP-1000	Cryo-Dipper	1	£118
CNET-1000	Cryo-Net	5	£24

In association with... **DIVERSIFIED BIOTECH**