


Quick-Therm

Thermally Conductive Aluminium Beads

Effectively heat or cool your samples in most water/dry baths, ice pans/buckets using our revolutionary Quick-Therm Conductive Aluminium Beads. Like our Thermal Microbeads (Prod Code: BEAD-2000) Quick-Therm beads are sterilizable and re-usable making them a convenient addition to any busy laboratory.



FEATURES

- Rapidly cool or heat samples in either plastic or glass vessels
- Aluminium mini-pellets thermally equilibrate samples up to 40% faster than larger beads.
- Samples are kept free of contaminants by eliminating water and ice
- Pellets conveniently support test tubes, microtubes, beakers and flasks without racks, floats or weights
- Pellets can be cleaned and disinfected with a 70% alcohol rinse

Order Information

Cat Number	Description	Pack Size	List Price
BEAD-1000	Quick-Therm Mini-Pellets	1 Litre	£47